

Face Scan Report

Stroke Detection Screening Analysis by NeuroScan

Scan ID: NS-202606150852
Scan Date: June 15, 2026
Scan Time: 08:52 AM
Methodology: MediaPipe Face Mesh + ML Classifier


OVERALL ASYMMETRY SCORE

99.6%

HIGH RISK — SEEK IMMEDIATE
MEDICAL HELP

Visual Facial Mesh Examination

Below is the client's photograph annotated with the MediaPipe 468-point Face Mesh. The system has highlighted landmarks where facial asymmetry was analyzed. Green circles indicate symmetric regions within normal boundaries, while red circles mark landmarks that exceed clinical asymmetry thresholds.




Detailed Asymmetry Breakdown

The table below details the 6 key facial asymmetry metrics extracted by MediaPipe landmarks. Each score represents the normalized difference between the left and right sides of the face (0 to 1 scale). A feature is flagged as **HIGH** if its measured value exceeds the clinically established threshold.

Facial Feature Metric	Measured Value	Normal Threshold	Clinical Basis / Reference	Symmetry Status
Eye Asymmetry (Ptosis)	0.281	≤ 0.12	Ross et al. (Sunnybrook FGS)	✗ HIGH
Mouth Corner Droop	0.043	≤ 0.10	Kothari et al. (CPSS FAST)	✓ Normal
Eyebrow Height Difference	0.363	≤ 0.08	Ross et al. (Sunnybrook FGS)	✗ HIGH
Nasolabial Fold Asymmetry	0.101	≤ 0.10	Taufique & Savakis (2021)	✗ HIGH
Facial Midline Deviation	0.062	≤ 0.05	Taufique & Savakis (2021)	✗ HIGH
Mouth Width Asymmetry	0.000	≤ 0.10	Kothari et al. (CPSS FAST)	✓ Normal

Risk Level Calibration



CRITICAL FINDING: Significant facial asymmetry has been detected in key facial regions, exceeding established clinical thresholds. In accordance with the FAST protocol, this is a strong indicator of potential acute stroke or facial nerve palsy.

RECOMMENDED ACTION: Please seek emergency medical help immediately. Call emergency services (**999 or 911**) or go to the nearest emergency department without delay. Do not wait to see if symptoms improve.

Clinical References

1. Herpich et al. Human vs. Machine Learning Based Detection of Facial Weakness Using Video Analysis. *Front. Neurol.*, 2022.
2. Taufique, A. & Savakis, A. Automatic Quantification of Facial Asymmetry Using Facial Landmarks. *arXiv:2103.11059*, 2021.
3. Kothari et al. Cincinnati Prehospital Stroke Scale: reproducibility and validity. *Acad. Emerg. Med.*, 1999.

Medical Disclaimer: This report is generated automatically by an artificial intelligence screening tool. It is designed solely for screening purposes to identify potential facial asymmetry. This report does NOT constitute medical advice, a diagnosis, or a treatment plan. A negative scan does not guarantee the absence of a stroke, and a positive scan does not guarantee the presence of a stroke. Always consult with a licensed healthcare professional for any clinical evaluation or diagnosis.